



EURL- *Salmonella*

Proficiency test on the detection of *Salmonella* in boot socks with chicken faeces (2018)



Irene Pol-Hofstad

EURL- Salmonella | 29 May 2019

EURL-*Salmonella* Detection study PPS 2018

- Matrix: Boot socks with chicken faeces
- Boot socks are used to sample poultry farms
- Boot socks with chicken faeces contaminated at EURL
 - high and low concentration of *Salmonella* Infantis




Preparation of samples

- 10 ml Peptone Saline Solution to moisten boot socks
- 15 g of chicken faeces
- Add high and low concentration *Salmonella* Infantis
 - Low: approx. 5-10 cfu
 - High: approx. 50-100 cfu



Pre-test: Stability of samples

- Simulate transport and storage conditions
- Storage conditions: 5 °C and 10 °C for up to 3 weeks
- Test for:
 - *Salmonella*: EN-ISO 6579:1-2017
 - *Enterobacteriaceae*: EN-ISO 21528-2
 - Aerobic count: EN-ISO 4833-1:2013




Conclusion Pre-tests: *Salmonella*

– *Salmonella*

> # positive samples


14 cfu



Dark: 5°C Light: 10°C

Conclusion Pre-tests: background flora

- *Aerobic count*
- *Enterobacteriaceae*



Study design

- Boot sock + 10 ml PSS
- 15 g chicken faeces
- Contaminated with *Salmonella* Infantis:
 - 6 x Low level SI: 14 cfu
 - 6 x High level SI: 50 cfu
 - 6 x Negative samples (No *Salmonella* added)
- Process controls:
 - C1: Boot sock + 10 ml PSS
 - C2: Boot sock + 10 ml PSS+ own positive control





Participants

In total 36 participating NRLs

- 29 NRLs from 28 EU Member States
- 6 NRLs from EU (potential) candidate Member States or European Free Trade Associations (EFTA)
- 1 NRL from a third country





Transport


- Transport time:
 - 26 parcels: 0-1 day
 - 6 parcels: 2 days
 - 2 parcels: 3 days
 - 2 parcels: problems at customs
 - labcode 35: arrival after 8 days
 - labcode 22: arrival after 10 days

- Temperature during transport and storage:
 - temperature probe between samples
 - transport temperature: predominantly $-3\text{ }^{\circ}\text{C}$ and $+5\text{ }^{\circ}\text{C}$
 - storage temperature: mostly $0\text{ }^{\circ}\text{C}$ - $6\text{ }^{\circ}\text{C}$.






Temperature samples parcel Lab 35



Method PT-PPS

Prescribed: ISO 6579-1:2017

- Selective enrichment on MSRV
- Isolation medium XLD & 2nd agar of choice





PT samples: Background flora

| Date of testing | Aerobic bacteria cfu/g | <i>Enterobacteriaceae</i> cfu/g |
|---|---------------------------|------------------------------------|
| Preparation: t = 0 days 17 Sept 2018 | 4.9×10^8 | 2.8×10^6 |
| Start analysis t = 14 days 1 Oct 2018 (after storage at 5 °C) | 1.2×10^8 | 1.9×10^5 |

PT samples: Artificial contamination

| Date of testing | Sample | Low level SI cfu/sample | High level SI cfu/sample |
|-----------------|---|----------------------------|-----------------------------|
| 18 Sept 2018 | Inoculum of diluted culture | 10 | 53 |
| 1 Oct 2018 | Inoculated boot socks with chicken faeces, stored 5°C MPN (95% confidence limit) | 3.3 (1.1 - 10.3) | 17.3 (6.5 - 45) |




Results: High concentration SI





Results Low concentration SI





Samples: specificity, sensitivity, accuracy

| Boot sock samples with chicken faeces | | All labs n = 35 | EU-MS n = 29 |
|--|-------------------------|--------------------|-----------------|
| Negative samples n=6 | No. of samples | 210 | 174 |
| | No. of negative samples | 210 | 174 |
| | Specificity in % | 100% | 100% |
| Low level (SI) n=6 | No. of samples | 210 | 174 |
| | No. of positive samples | 208 | 173 |
| | Sensitivity in % | 99% | 99.4% |
| High level (SI) n=6 | No. of samples | 210 | 174 |
| | No. of positive samples | 209 | 173 |
| | Sensitivity in % | 99.5% | 99.4% |
| All boot sock samples with SI | No. of samples | 420 | 348 |
| | No. of positive samples | 417 | 346 |
| | Sensitivity in % | 99,3% | 99.4% |
| All boot sock samples (positive and negative) | No. of samples | 630 | 522 |
| | No. of correct samples | 627 | 520 |
| | Accuracy in % | 99,5% | 99,6% |



Control samples: specificity, sensitivity, accuracy

| Control samples | | All labs n = 36 | EU n = 29 |
|---|-------------------------|--------------------|--------------|
| Procedure control (moistened boot socks) n = 1 | No. of samples | 36 | 29 |
| | No. of negative samples | 36 | 29 |
| | Specificity in % | 100% | 100% |
| Positive control (Own <i>Salmonella</i>) n = 1 | No. of samples | 36 | 29 |
| | No. of positive samples | 36 | 29 |
| | Sensitivity in % | 100% | 100% |
| All control samples n = 2 | No. of samples | 72 | 58 |
| | No. of correct samples | 72 | 58 |
| | Accuracy in % | 100% | 100% |



Conclusions PT-PPS 2018

- Control samples:
 - Excellent performance: 100% correct scores
- *Salmonella* boot sock samples with chicken faeces:
 - High level SI: 1 lab scored 1 sample negative
 - Low level SI: 2 labs scores 1 sample negative
 - Negative samples: all labs scored all samples negative
- Results of lab 35 were not evaluated because of temperature abuse during transport



Thank you all for your participation in this study

