


**Example of the on-line Result form
EURL-Salmonella Proficiency Test PPS 2019**

k = Selecteren
klk = Menu
clk = Bewerken



EURL-*Salmonella* Proficiency Test PPS 2019

Detection of *Salmonella* in chicken faeces

Result Form

Laboratory information

Laboratory code

Name contact person

E-mail address contact person

Name laboratory or institute

Country

For which *Salmonella* method(s) is your laboratory accredited?

EN ISO 6579:2002

Annex D of EN ISO 6579:2007

EN ISO 6579-1:2017

Other. Please specify:

Which method did you use in this PT?

EN ISO 6579:2002

Annex D of EN ISO 6579:2007

EN ISO 6579-1:2017

Other name:

Date and time of arrival of the parcel in your laboratory

dd/mm/yyyy Time

Was your parcel damaged at arrival?

No
 Yes, Description of damage:

Start date and time of storage at 5°C

dd/mm/yyyy Time of storage

Start date of proficiency test

dd/mm/yyyy

Positive control sample

Which *Salmonella* serovar did you use in the control samples?

- Enteritidis
- Typhimurium
- Goldcoast
- Nottingham
- Panama
- Poona
- Tranaroa (*Salmonella enterica* subsp. *salamae*)
- Other

What was the concentration of the *Salmonella* control sample?

Was matrix added to the positive control sample?

- No
- Yes, specify the matrix:

Pre-enrichment - Buffered Peptone Water (BPW)

Temperature incubator BPW

°C

Hours of incubation

Selective enrichment - Modified Semi solid Rappaport medium (MSRV)

Concentration
novobiocin per 1L
medium

Temperature
incubation MSRV

pH at the day of
use

Selective Isolation media

Which isolation
media were
used?

- ASAP (AES Salmonella Agar Plate)
- BGA
- BGA (mod)
- BSA
- BPLS
- BxLH
- Rambach
- SM(ID)2
- XLD
- Other, please specify

Confirmation of *Salmonella* suspected colonies

What type of tests
did you use for
confirmation?

- Biochemical
- Serological
- Serotyping
- PCR
- Other Specify:

Second Detection Method (optional)

Were the samples also tested with a second method? Yes No

What type of method? For example: PCR, qPCR, Rapid *Salmonella*

Is this method validated? Yes No

Which organisation validated the method?

Please give a certificate number of reference for the validated method

Is the second method used for analysing official control samples? Yes No

How many samples did you (approx.) analyse with this second method in 2018?

Table 1: Confirmed *Salmonella* results

Indicate:

detected for positive confirmed result for *Salmonella*

not detected for negative (confirmed) results for *Salmonella*

Method → ISO 6579-1

Second method

Samples		
↓	B1	B1.
	B2	B2.
	B3	B3.
	B4	B4.
	B5	B5.
	B6	B6.
	B7	B7.
	B8	B8.
	B9	B9.
	B10	B10.
	B11	B11.
	B12	B12.
	B13	B13.
	B14	B14.
	C1	C1.
	C2	C2.

Remarks and comments

Name of person(s) carrying out the comparison study

Name of the person in charge

SAVE YOUR PROGRESS RESUME LATER

TO CONFIRMATION PAGE >>